





CYP1B1 Antibody

Product Code	CSB-PA001931
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16678
Immunogen	Synthesized peptide derived from the Internal region of Human CYP1B1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CYP1B1; Cytochrome P450 1B1; CYPIB1
Immunogen Species	Homo sapiens (Human)
Target Names	CYP1B1

Image



Western Blot analysis of Jurkat cells using CYP1B1 Polyclonal Antibody



Western Blot analysis of Jurkat cells using CYP1B1 Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.