





CYP2A7 Antibody

Product Code	CSB-PA001953
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20853
Immunogen	Synthesized peptide derived from the C-terminal region of Human CYP2A7.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CYP2A7; Cytochrome P450 2A7; CYP11A7; Cytochrome P450 IIA4
Immunogen Species	Homo sapiens (Human)
Target Names	CYP2A7

Image



Western Blot analysis of 22RV1 cells using CYP2A7 Polyclonal Antibody



Western Blot analysis of Jurkat cells using CYP2A7 Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.