




CYP2B6 Antibody

Product Code	CSB-PA001956
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20813
Immunogen	Synthesized peptide derived from the Internal region of Human CYP2B6.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CYP2B6; Cytochrome P450 2B6; 1; 4-cineole 2-exo-monooxygenase; CYP11B6; Cytochrome P450 IIB1
Immunogen Species	Homo sapiens (Human)
Target Names	CYP2B6

Image



Western Blot analysis of HT29 cells using CYP2B6 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.