




CYP2W1 Antibody

Product Code	CSB-PA001971
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAV3
Immunogen	Synthesized peptide derived from the Internal region of Human CYP2W1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CYP2W1; Cytochrome P450 2W1; CYP1IW1
Immunogen Species	Homo sapiens (Human)
Target Names	CYP2W1

Image



Western Blot analysis of HepG2 cells using CYP2W1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.