





CYP39A1 Antibody

Product Code	CSB-PA001977
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYL5
Immunogen	Synthesized peptide derived from the Internal region of Human CYP39A1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CYP39A1; 24-hydroxycholesterol 7-alpha-hydroxylase; Cytochrome P450 39A1; hCYP39A1; Oxysterol 7-alpha-hydroxylase
Immunogen Species	Homo sapiens (Human)
Target Names	CYP39A1

Image



Western Blot analysis of HeLa cells using CYP39A1 Polyclonal Antibody



Western Blot analysis of K562 cells using
CYP39A1 Polyclonal Antibody