




CYP4Z2P Antibody

Product Code	CSB-PA001987
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VVE4/Q8N1L4
Immunogen	Synthesized peptide derived from the Internal region of Human CYP4Z1/2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	Cytochrome P450 4Z1/Putative cytochrome P450 family member 4Z2
Immunogen Species	Homo sapiens (Human)
Target Names	CYP4Z2P

Image



Western Blot analysis of Jurkat RAT-MUSCLE cells using CYP4Z1/2 Polyclonal Antibody