




CST11 Antibody

Product Code	CSB-PA001994
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H112
Immunogen	Synthesized peptide derived from the C-terminal region of Human Cystatin 11.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CST11; CST8L; Cystatin-11
Immunogen Species	Homo sapiens (Human)
Target Names	CST11

Image



Western Blot analysis of A549 cells using Cystatin 11 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.