





CST3 Antibody

Product Code	CSB-PA002010
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01034
Immunogen	Synthesized peptide derived from the Internal region of Human Cystatin C.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CST3; Cystatin-C; Cystatin-3; Gamma-trace; Neuroendocrine basic polypeptide; Post-gamma-globulin
Immunogen Species	Homo sapiens (Human)
Target Names	CST3

Image



Western Blot analysis of 293T cells using Cystatin C Polyclonal Antibody



Western Blot analysis of 293T cells using Cystatin C Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.