





KRT18 Antibody

Product Code	CSB-PA002029
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05783
Immunogen	Synthesized peptide derived from Human Cytokeratin 18 around the non-phosphorylation site of S33.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IP, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	KRT18; CYK18; PIG46; Keratin; type I cytoskeletal 18; Cell proliferation-inducing gene 46 protein; Cytokeratin-18; CK-18; Keratin-18; K18
Immunogen Species	Homo sapiens (Human)
Target Names	KRT18

Image



Western Blot analysis of HT29 cells using Cytokeratin 18 Polyclonal Antibody



Western Blot analysis of VEC cells using Cytokeratin 18 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.