




KRT20 Antibody

Product Code	CSB-PA002048
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35900
Immunogen	Synthesized peptide derived from the C-terminal region of Human Cytokeratin 20.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	KRT20; Keratin; type I cytoskeletal 20; Cytokeratin-20; CK-20; Keratin-20; K20; Protein IT
Immunogen Species	Homo sapiens (Human)
Target Names	KRT20

Image



Western Blot analysis of HeLa cells using Cytokeratin 20 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.