




TPD52L2 Antibody

Product Code	CSB-PA002059
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43399
Immunogen	Synthesized peptide derived from the C-terminal region of Human D54.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TPD52L2; Tumor protein D54; hD54; Tumor protein D52-like 2
Immunogen Species	Homo sapiens (Human)
Target Names	TPD52L2

Image



Western Blot analysis of Jurkat cells using D54 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.