





DAB1 Antibody

Product Code	CSB-PA002060
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75553
Immunogen	Synthesized peptide derived from Human Dab1 around the non-phosphorylation site of Y232.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DAB1; Disabled homolog 1
Immunogen Species	Homo sapiens (Human)
Target Names	DAB1

Image



Western Blot analysis of 293 cells using Dab1 Polyclonal Antibody



Western Blot analysis of RAT-MUSCLE cells using Dab1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.