




NBL1 Antibody

Product Code	CSB-PA002066
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41271
Immunogen	Synthesized peptide derived from the C-terminal region of Human DAN.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NBL1; DAN; DAND1; Neuroblastoma suppressor of tumorigenicity 1; DAN domain family member 1; Protein N03; Zinc finger protein DAN
Immunogen Species	Homo sapiens (Human)
Target Names	NBL1

Image



Western Blot analysis of COLO205 cells using DAN Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.