





PPP1R1B Antibody

Product Code	CSB-PA002075
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UD71
Immunogen	Synthesized peptide derived from Human DARPP-32 around the non-phosphorylation site of T75.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PPP1R1B; DARPP32; Protein phosphatase 1 regulatory subunit 1B; DARPP-32; Dopamine- and cAMP-regulated neuronal phosphoprotein
Immunogen Species	Homo sapiens (Human)
Target Names	PPP1R1B

Image



Western Blot analysis of COLO205 cells using DARPP-32 Polyclonal Antibody



Western Blot analysis of various cells using DARPP-32 Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.