





DAXX Antibody

Product Code	CSB-PA002077
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UER7
Immunogen	Synthesized peptide derived from the Internal region of Human Daxx.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DAXX; BING2; DAP6; Death domain-associated protein 6; Daxx; hDaxx; ETS1-associated protein 1; EAP1; Fas death domain-associated protein
Immunogen Species	Homo sapiens (Human)
Target Names	DAXX

Image



Western Blot analysis of HuvEc cells using Daxx Polyclonal Antibody



Western Blot analysis of 293 cells using Daxx Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.