




DDX51 Antibody

Product Code	CSB-PA002107
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N8A6
Immunogen	Synthesized peptide derived from the C-terminal region of Human DDX51.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DDX51; ATP-dependent RNA helicase DDX51; DEAD box protein 51
Immunogen Species	Homo sapiens (Human)
Target Names	DDX51

Image



Western Blot analysis of HepG2 cells using DDX51 Polyclonal Antibody