




DDX52 Antibody

Product Code	CSB-PA002108
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2R4
Immunogen	Synthesized peptide derived from the C-terminal region of Human DDX52.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DDX52; ROK1; HUSSY-19; Probable ATP-dependent RNA helicase DDX52; ATP-dependent RNA helicase ROK1-like; DEAD box protein 52
Immunogen Species	Homo sapiens (Human)
Target Names	DDX52

Image



Western Blot analysis of HepG2 cells using DDX52 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.