




DDX55 Antibody

Product Code	CSB-PA002112
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NHQ9
Immunogen	Synthesized peptide derived from the Internal region of Human DDX55.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DDX55; KIAA1595; ATP-dependent RNA helicase DDX55; DEAD box protein 55
Immunogen Species	Homo sapiens (Human)
Target Names	DDX55

Image



Western Blot analysis of HeLa cells using DDX55 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.