





DGKA Antibody

Product Code	CSB-PA002124
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23743
Immunogen	Synthesized peptide derived from the Internal region of Human DGK- α .
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DGKA; DAGK; DAGK1; Diacylglycerol kinase alpha; DAG kinase alpha; 80 kDa diacylglycerol kinase; Diglyceride kinase alpha; DGK-alpha
Immunogen Species	Homo sapiens (Human)
Target Names	DGKA

Image



Western Blot analysis of Lovo cells using DGK- α Polyclonal Antibody



Western Blot analysis of LOVO cells using DGK- α Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.