




DGKD Antibody

Product Code	CSB-PA002129
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16760
Immunogen	Synthesized peptide derived from the N-terminal region of Human DGK- δ .
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DGKD; KIAA0145; Diacylglycerol kinase delta; DAG kinase delta; 130 kDa diacylglycerol kinase; Diglyceride kinase delta; DGK-delta
Immunogen Species	Homo sapiens (Human)
Target Names	DGKD

Image



Western Blot analysis of HuvEc cells using DGK- δ Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.