





DHPS Antibody

Product Code	CSB-PA002145
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49366
Immunogen	Synthesized peptide derived from the Internal region of Human DHS.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DHPS; DS; Deoxyhypusine synthase; DHS
Immunogen Species	Homo sapiens (Human)
Target Names	DHPS

Image



Western Blot analysis of K562 cells using DHS Polyclonal Antibody



Western Blot analysis of K562 cells using DHS Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.