





DMPK Antibody

Product Code	CSB-PA002164
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q09013
Immunogen	Synthesized peptide derived from the N-terminal region of Human DMPK.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DMPK; DM1PK; MDPK; Myotonin-protein kinase; MT-PK; DM-kinase; DMK; DM1 protein kinase; DMPK; Myotonic dystrophy protein kinase
Immunogen Species	Homo sapiens (Human)
Target Names	DMPK

Image



Western Blot analysis of Jurkat cells using DMPK Polyclonal Antibody



Western Blot analysis of Jurkat cells using DMPK Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.