




LIG1 Antibody

Product Code	CSB-PA002165
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18858
Immunogen	Synthesized peptide derived from the N-terminal region of Human DNA Ligase I.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	LIG1; DNA ligase 1; DNA ligase I; Polydeoxyribonucleotide synthase [ATP] 1
Immunogen Species	Homo sapiens (Human)
Target Names	LIG1

Image



Western Blot analysis of HT29 cells using DNA Ligase I Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.