





LIG4 Antibody

Product Code	CSB-PA002168
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49917
Immunogen	Synthesized peptide derived from the Internal region of Human DNA Ligase IV.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	LIG4; DNA ligase 4; DNA ligase IV; Polydeoxyribonucleotide synthase [ATP] 4
Immunogen Species	Homo sapiens (Human)
Target Names	LIG4

Image



Western Blot analysis of VEC cells using DNA Ligase IV Polyclonal Antibody



Western Blot analysis of 293 A549 cells using
DNA Ligase IV Polyclonal Antibody