





LIG4 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA002168 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P49917 |
| Immunogen | Synthesized peptide derived from the Internal region of Human DNA Ligase IV. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | LIG4; DNA ligase 4; DNA ligase IV; Polydeoxyribonucleotide synthase [ATP] 4 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | LIG4 |

Image



Western Blot analysis of VEC cells using DNA Ligase IV Polyclonal Antibody



Western Blot analysis of 293 A549 cells using DNA Ligase IV Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.