




POLL Antibody

Product Code	CSB-PA002174
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGP5
Immunogen	Synthesized peptide derived from the C-terminal region of Human DNA pol λ.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	POLL; DNA polymerase lambda; Pol Lambda; DNA polymerase beta-2; Pol beta2; DNA polymerase kappa
Immunogen Species	Homo sapiens (Human)
Target Names	POLL

Image



Western Blot analysis of 293 cells using POLL Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.