




DNMT3B Antibody

Product Code	CSB-PA002178
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBC3
Immunogen	Synthesized peptide derived from the N-terminal region of Human Dnmt3b.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DNMT3B; DNA (cytosine-5)-methyltransferase 3B; Dnmt3b; DNA methyltransferase HsallIB; DNA MTase HsallIB; M.HsallIB
Immunogen Species	Homo sapiens (Human)
Target Names	DNMT3B

Image



Western Blot analysis of HT29 cells using Dnmt3b Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.