




TFDP2 Antibody

Product Code	CSB-PA002187
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14188
Immunogen	Synthesized peptide derived from the Internal region of Human DP-2.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TFDP2; DP2; Transcription factor Dp-2; E2F dimerization partner 2
Immunogen Species	Homo sapiens (Human)
Target Names	TFDP2

Image



Western Blot analysis of COS7 cells using DP-2 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.