




DQX1 Antibody

Product Code	CSB-PA002190
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TE96
Immunogen	Synthesized peptide derived from the C-terminal region of Human DQX1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DQX1; ATP-dependent RNA helicase DQX1; DEAQ box polypeptide 1
Immunogen Species	Homo sapiens (Human)
Target Names	DQX1

Image



Western Blot analysis of COS7 cells using DQX1 Polyclonal Antibody