





STK17B Antibody

Product Code	CSB-PA002192
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94768
Immunogen	Synthesized peptide derived from the Internal region of Human DRAK2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	STK17B; DRAK2; Serine/threonine-protein kinase 17B; DAP kinase-related apoptosis-inducing protein kinase 2
Immunogen Species	Homo sapiens (Human)
Target Names	STK17B

Image



Western Blot analysis of brain cells using DRAK2 Polyclonal Antibody



Western Blot analysis of Jurkat K562 cells using DRAK2 Polyclonal Antibody