




GPR126 Antibody

Product Code	CSB-PA002202
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86SQ4
Immunogen	Synthesized peptide derived from the C-terminal region of Human DREG.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPR126; DREG; VIGR; G-protein coupled receptor 126; Developmentally regulated G-protein-coupled receptor; Vascular inducible G protein-coupled receptor
Immunogen Species	Homo sapiens (Human)
Target Names	ADGRG6

Image



Western Blot analysis of HepG2 cells using DREG Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.