





DYRK1A Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA002232 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q13627 |
| Immunogen | Synthesized peptide derived from the N-terminal region of Human Dyrk1A. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | DYRK1A; DYRK; MNB; MNBH; Dual specificity tyrosine-phosphorylation-regulated kinase 1A; Dual specificity YAK1-related kinase; HP86; Protein kinase minibrain homolog; MNBH; hMNB |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | DYRK1A |

Image



Western Blot analysis of HepG2 cells using Dyrk1A Polyclonal Antibody



Western Blot analysis of VEC cells using Dyrk1A Polyclonal Antibody