





DYRK1A Antibody

Product Code	CSB-PA002232
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13627
Immunogen	Synthesized peptide derived from the N-terminal region of Human Dyrk1A.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DYRK1A; DYRK; MNB; MNBH; Dual specificity tyrosine-phosphorylation-regulated kinase 1A; Dual specificity YAK1-related kinase; HP86; Protein kinase minibrain homolog; MNBH; hMNB
Immunogen Species	Homo sapiens (Human)
Target Names	DYRK1A

Image



Western Blot analysis of HepG2 cells using Dyrk1A Polyclonal Antibody



Western Blot analysis of VEC cells using Dyrk1A Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.