




E2F4/E2F5 Antibody

Product Code	CSB-PA002237
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16254/Q15329
Immunogen	Synthesized peptide derived from the Internal region of Human E2F-4/5.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	E2F4; Transcription factor E2F4; E2F-4; E2F5; Transcription factor E2F5; E2F-5
Immunogen Species	Homo sapiens (Human)
Target Names	E2F4/E2F5

Image



Western Blot analysis of COS7 cells using E2F-4/5 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.