




NR2F6 Antibody

Product Code	CSB-PA002242
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10588
Immunogen	Synthesized peptide derived from the N-terminal region of Human EAR2.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NR2F6; EAR2; ERBAL2; Nuclear receptor subfamily 2 group F member 6; V-erbA-related protein 2; EAR-2
Immunogen Species	Homo sapiens (Human)
Target Names	NR2F6

Image



Western Blot analysis of 293 cells using EAR2 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.