





CDH1 Antibody

Product Code	CSB-PA002243
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12830
Immunogen	Synthesized peptide derived from the N-terminal region of Human E-cadherin.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CDH1

Image



Western Blot analysis of 293 cells using E-cadherin Polyclonal Antibody



Western Blot analysis of COLO 293T HELA cells using E-cadherin Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.