




RNASE3 Antibody

Product Code	CSB-PA002246
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12724
Immunogen	Synthesized peptide derived from the C-terminal region of Human ECP.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RNASE3; ECP; RNS3; Eosinophil cationic protein; ECP; Ribonuclease 3; RNase 3
Immunogen Species	Homo sapiens (Human)
Target Names	RNASE3

Image



Western Blot analysis of 293 cells using ECP Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.