





LPAR3 Antibody

Product Code	CSB-PA002255
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBY5
Immunogen	Synthesized peptide derived from the C-terminal region of Human EDG-7.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	LPAR3; EDG7; LPA3; Lysophosphatidic acid receptor 3; LPA receptor 3; LPA-3; Lysophosphatidic acid receptor Edg-7
Immunogen Species	Homo sapiens (Human)
Target Names	LPAR3

Image



Western Blot analysis of 3T3 cells using EDG-7 Polyclonal Antibody



Western Blot analysis of various cells using EDG-7 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.