





EEF2K Antibody

Product Code	CSB-PA002258
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00418
Immunogen	Synthesized peptide derived from Human eEF2K around the non-phosphorylation site of S366.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EEF2K; Eukaryotic elongation factor 2 kinase; eEF-2 kinase; eEF-2K; Calcium/calmodulin-dependent eukaryotic elongation factor 2 kinase
Immunogen Species	Homo sapiens (Human)
Target Names	EEF2K

Image



Western Blot analysis of HeLa cells using eEF2K Polyclonal Antibody



Western Blot analysis of various cells using eEF2K Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.