




MEGF9 Antibody

Product Code	CSB-PA002263
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1U4
Immunogen	Synthesized peptide derived from the Internal region of Human EGFL5.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MEGF9; EGFL5; KIAA0818; Multiple epidermal growth factor-like domains protein 9; Multiple EGF-like domains protein 9; Epidermal growth factor-like protein 5; EGF-like protein 5
Immunogen Species	Homo sapiens (Human)
Target Names	MEGF9

Image



Western Blot analysis of COLO205 cells using EGFL5 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.