





EIF2S1 Antibody

Product Code	CSB-PA002295
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05198
Immunogen	Synthesized peptide derived from Human eIF2 α around the non-phosphorylation site of S51.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EIF2S1; EIF2A; Eukaryotic translation initiation factor 2 subunit 1; Eukaryotic translation initiation factor 2 subunit alpha; eIF-2-alpha; eIF-2A; eIF-2alpha
Immunogen Species	Homo sapiens (Human)
Target Names	EIF2S1

Image



Western Blot analysis of 3T3 cells using eIF2 α Polyclonal Antibody



Western Blot analysis of K562 cells using eIF2α Polyclonal Antibody