




EIF4E Antibody

Product Code	CSB-PA002308
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06730
Immunogen	Synthesized peptide derived from Human eIF4E around the non-phosphorylation site of S209.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EIF4E; EIF4EL1; EIF4F; Eukaryotic translation initiation factor 4E; eIF-4E; eIF4E; eIF-4F 25 kDa subunit; mRNA cap-binding protein
Immunogen Species	Homo sapiens (Human)
Target Names	EIF4E

Image



Western Blot analysis of COLO205 cells using eIF4E Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.