





KLF1/KLF5/KLF7 Antibody

Product Code	CSB-PA002315
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13351/Q13887/O75840
Immunogen	Synthesized peptide derived from the C-terminal region of Human EKLF/CKLF/UKLF.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	KLF1; EKLF; Krueppel-like factor 1; Erythroid krueppel-like transcription factor; EKLF; KLF5; BTEB2; CKLF; IKLF; Krueppel-like factor 5; Basic transcription element-binding protein 2; BTE-binding protein 2; Colon krueppel-like factor; GC-bo
Immunogen Species	Homo sapiens (Human)
Target Names	KLF1/KLF5/KLF7

Image



Western Blot analysis of VEC cells using EKLF/CKLF/UKLF Polyclonal Antibody



Western Blot analysis of VEC HY926 cells using EKLF/CKLF/UKLF Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.