




ELF5 Antibody

Product Code	CSB-PA002318
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKW6
Immunogen	Synthesized peptide derived from the C-terminal region of Human Elf-5.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ELF5; ESE2; ETS-related transcription factor Elf-5; E74-like factor 5; Epithelium-restricted ESE-1-related Ets factor; Epithelium-specific Ets transcription factor 2; ESE-2
Immunogen Species	Homo sapiens (Human)
Target Names	ELF5

Image



Western Blot analysis of HuvEc cells using Elf-5 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.