




ELOVL4 Antibody

Product Code	CSB-PA002325
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZR5
Immunogen	Synthesized peptide derived from the Internal region of Human ELOVL4.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ELOVL4; Elongation of very long chain fatty acids protein 4; 3-keto acyl-CoA synthase ELOVL4; ELOVL fatty acid elongase 4; ELOVL FA elongase 4
Immunogen Species	Homo sapiens (Human)
Target Names	ELOVL4

Image



Western Blot analysis of HeLa cells using ELOVL4 Polyclonal Antibody