




EOMES Antibody

Product Code	CSB-PA002351
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95936
Immunogen	Synthesized peptide derived from the N-terminal region of Human EOMES.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EOMES; TBR2; Eomesodermin homolog; T-box brain protein 2; T-brain-2; TBR-2
Immunogen Species	Homo sapiens (Human)
Target Names	EOMES

Image



Western Blot analysis of HepG2 cells using EOMES Polyclonal Antibody