




EOMES Antibody

Product Code	CSB-PA002351
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95936
Immunogen	Synthesized peptide derived from the N-terminal region of Human EOMES.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EOMES; TBR2; Eomesodermin homolog; T-box brain protein 2; T-brain-2; TBR-2
Immunogen Species	Homo sapiens (Human)
Target Names	EOMES

Image



Western Blot analysis of HepG2 cells using EOMES Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.