





PTGER4 Antibody

Product Code	CSB-PA002354
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35408
Immunogen	Synthesized peptide derived from the Internal region of Human EP4.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PTGER4; PTGER2; Prostaglandin E2 receptor EP4 subtype; PGE receptor EP4 subtype; PGE2 receptor EP4 subtype; Prostanoid EP4 receptor
Immunogen Species	Homo sapiens (Human)
Target Names	PTGER4

Image



Western Blot analysis of COLO205 cells using EP4 Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using EP4 Polyclonal Antibody