




EFNB1 Antibody

Product Code	CSB-PA002376
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98172
Immunogen	Synthesized peptide derived from Human Ephrin-B1 around the non-phosphorylation site of Y317.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EFNB1; EFL3; EPLG2; LERK2; Ephrin-B1; EFL-3; ELK ligand; ELK-L; EPH-related receptor tyrosine kinase ligand 2; LERK-2
Immunogen Species	Homo sapiens (Human)
Target Names	EFNB1

Image



Western Blot analysis of COS7 cells using Ephrin-B1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.