





ETV1 Antibody

Product Code	CSB-PA002391
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50549
Immunogen	Synthesized peptide derived from the N-terminal region of Human ER81.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ETV1; ER81; ETS translocation variant 1; Ets-related protein 81
Immunogen Species	Homo sapiens (Human)
Target Names	ETV1

Image



Western Blot analysis of HuvEc cells using ER81 Polyclonal Antibody



Western Blot analysis of K562 RAT-intestine cells using ER81 Polyclonal Antibody

